

Product Datasheet

IFT172 antibody (FITC) (orb516237)

Description	IFT172 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Intraflagellar transport protein 172 homolog protein (1368-1502AA)
Target	IFT172
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9UG01
Expiration Date	12 months from date of receipt.