

Product Datasheet

2-hydroxyisobutyryl-HIST1H4A (K16) antibody (orb516315)

Description 2-hydroxyisobutyryl-HIST1H4A (K16) antibody

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications ChIP, ELISA, IF, IP, WB

Immunogen Peptide sequence around site of 2-hydroxyisobutyryl-Lys (16) derived from Human Histone H4

Target HIST1H4A

Preservatives Preservative: 0.03% Proclin 300.
Constituents: 50% Glycerol

Form/Appearance Liquid

Storage Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.

Note For research use only


Isotype IgG

Clonality Polyclonal


Uniprot ID **P62805**

Dilution Range WB:1:500-1:5000, IF:1:50-1:200, IP:1:200-1:2000


Expiration Date 12 months from date of receipt.



Chromatin Immunoprecipitation
Hela (4*10...



Western blot analysis of Hela whole cel...



Immunoprecipitation analysis of HepG2 wh...