



Product Datasheet

2-hydroxyisobutyryl-HIST1H4A (K5) antibody (orb516320)


Description	2-hydroxyisobutyryl-HIST1H4A (K5) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, WB
Immunogen	Peptide sequence around site of 2-hydroxyisobutyryl-Lys (5) derived from Human Histone H4
Target	HIST1H4A
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P62805
Dilution Range	WB:1:500-1:5000, ICC:1:20-1:200, IF:1:50-1:200
Expiration Date	12 months from date of receipt.



Immunocytochemistry analysis of HeLa cells...



Western blot analysis of 293 whole cell...



Immunofluorescence staining of HeLa cells...