



Product Datasheet

HIST1H4A (Ab-77) antibody (orb516321)

Description	HIST1H4A (Ab-77) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Peptide sequence around site of Lys (77) derived from Human Histone H4
Target	HIST1H4A
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P62805
Dilution Range	WB:1:500-1:5000, IHC-P:1:1-1:10
Expiration Date	12 months from date of receipt.



Immunohistochemistry of paraffin-embedde...



Western blot aanalysis of HeLa whole cel...