



Product Datasheet

Crotonyl-HIST1H2BC (K23) antibody (orb516339)


Description	Crotonyl-HIST1H2BC (K23) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ChIP, ELISA, IF, IP, WB
Immunogen	Peptide sequence around site of Crotonyl-Lys (23) derived from Human Histone H2B type 1-C/E/F/G/I
Target	HIST1H2BC
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P62807
Dilution Range	WB:1:500-1:5000, IF:1:50-1:200, IP:1:200-1:2000
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of HeLa cell...



Western blot analysis of HeLa whole cel...



Immunoprecipitation analysis of K562 who...