



Product Datasheet

HIST1H4A (Ab-12) antibody (orb516346)


Description	HIST1H4A (Ab-12) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, IP, WB
Immunogen	Peptide sequence around site of Lys (12) derived from Human Histone H4
Target	HIST1H4A
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P62805
Dilution Range	WB:1:500-1:5000, IHC-P:1:20-1:200, IF:1:50-1:200, IP:1:200-1:2000
Expiration Date	12 months from date of receipt.



Immunohistochemistry of paraffin-embedde...



Western blot aanalysis of HeLa whole cel...



Immunoprecipitation analysis of HeLa who...