


Product Datasheet

Di-methyl-HIST1H2AG (R29) antibody (orb516385)

Description	Di-methyl-HIST1H2AG (R29) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ChIP, ELISA
Immunogen	Peptide sequence around site of Di-methyl-Arg (29) derived from Human Histone H2A type 1
Target	HIST1H2AG
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P0C0S8
Expiration Date	12 months from date of receipt.



Chromatin
Immunoprecipitation
Hela (4*10...