

Product Datasheet

KRIT1 antibody (FITC) (orb516460)

Description	KRIT1 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Krev interaction trapped protein 1 protein (471-585AA)
Target	KRIT1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	000522
Expiration Date	12 months from date of receipt.