

Product Datasheet

vjbR antibody (FITC) (orb51652)

Description	vjbR antibody (FITC)
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	FITC
Immunogen	Recombinant Brucella abortus HTH-type quorum sensing-dependent transcriptional regulator VjbR protein (1-259AA)
Target	vjbR
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q2YJ50
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.