

Product Datasheet

Guanine nucleotide-binding protein subunit beta-like protein antibody (Biotin) (orb51657)

Description	Guanine nucleotide-binding protein subunit beta-like protein antibody (Biotin)
Species/Host	Rabbit
Reactivity	Plant
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Glycine max Guanine nucleotide-binding protein subunit beta-like protein (1-325AA)
Target	Guanine nucleotide-binding protein subunit beta-like protein
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q39836
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.