

Product Datasheet

FREM1 antibody (FITC) (orb516679)

Description	FREM1 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human FRAS1-related extracellular matrix protein 1 protein (1957-2082AA)
Target	FREM1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q5H8C1
Expiration Date	12 months from date of receipt.