

Product Datasheet

ALG10 antibody (Biotin) (orb516686)

Description	ALG10 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human DoI-P-Glc:Glc(2)Man(9)GlcNAc(2)-PP-DoI alpha-1,2-glucosyltransferase protein (197-249AA)
Target	ALG10
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q5BKT4
Expiration Date	12 months from date of receipt.