

Product Datasheet

SLC6A19 antibody (Biotin) (orb516710)

Description	SLC6A19 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Sodium-dependent neutral amino acid transporter B(0)AT1 protein (326-413AA)
Target	SLC6A19
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q695T7
Expiration Date	12 months from date of receipt.