

Product Datasheet

NKX2-3 antibody (FITC) (orb516850)

Description	NKX2-3 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Homeobox protein Nkx-2.3 protein (52-148AA)
Target	NKX2-3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8TAU0
Expiration Date	12 months from date of receipt.