

Product Datasheet

MFSD4A antibody (FITC) (orb516853)

Description	MFSD4A antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Major facilitator superfamily domain-containing protein 4A protein (160-220AA)
Target	MFSD4A
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8N468
Expiration Date	12 months from date of receipt.