

Product Datasheet

INPP5E antibody (FITC) (orb516901)

Description	INPP5E antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human 72 kDa inositol polyphosphate 5-phosphatase protein (1-133AA)
Target	INPP5E
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9NRR6
Expiration Date	12 months from date of receipt.