

Product Datasheet

CXorf36 antibody (FITC) (orb516946)

Description	CXorf36 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Deleted in autism-related protein 1 protein (311-405AA)
Target	CXorf36
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9H7Y0
Expiration Date	12 months from date of receipt.