



Product Datasheet

Crotonyl-HIST1H2BC (K16) antibody (orb517017)


Description	Crotonyl-HIST1H2BC (K16) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ChIP, ELISA, ICC, IF, WB
Immunogen	Peptide sequence around site of Crotonyl-Lys (16) derived from Human Histone H2B type 1-C/E/F/G/I
Target	HIST1H2BC
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P62807
Dilution Range	WB:1:500-1:5000, ICC:1:500-1:1000, IF:1:50-1:200
Expiration Date	12 months from date of receipt.



Chromatin Immunoprecipitation
Hela (4*10...



Western blot
analysis of Hela
whole cel...



Immunofluorescence
staining of Hela
cell...