



Product Datasheet

2-hydroxyisobutyryl-HIST1H2BC (K12) antibody (orb517035)


Description	2-hydroxyisobutyryl-HIST1H2BC (K12) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, WB
Immunogen	Peptide sequence around site of 2-hydroxyisobutyryl-Lys (12) derived from Human Histone H2B type 1-C/E/F/G/I
Target	HIST1H2BC
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P62807
Dilution Range	WB:1:500-1:5000, ICC:1:200-1:500, IF:1:50-1:200
Expiration Date	12 months from date of receipt.



Immunocytochemistry analysis of Hela cel...



Western blot aanalysis of K562 whole cel...



Immunofluorescence staining of Hela cell...