

Product Datasheet

hlgC antibody (Biotin) (orb51705)

Description	hlgC antibody (Biotin)
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Staphylococcus aureus Gamma-hemolysin component C protein (30-315AA)
Target	hlgC
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q6GE13
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.