



Product Datasheet

Butyrlly-HIST1H3A (K18) antibody (orb517062)


Description	Butyrlly-HIST1H3A (K18) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, WB
Immunogen	Peptide sequence around site of Butyrlly-Lys (18) derived from Human Histone H3.1
Target	HIST1H3A
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P68431
Dilution Range	WB:1:500-1:5000, ICC:1:200-1:500, IF:1:50-1:200
Expiration Date	12 months from date of receipt.



Immunocytochemistry analysis of HeLa cell...



Western blot analysis of HeLa whole cel...



Immunofluorescence staining of HeLa cell...