


Product Datasheet

2-hydroxyisobutyryl-HIST1H2AG (K36) antibody (orb517081)

Description	2-hydroxyisobutyryl-HIST1H2AG (K36) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, ICC
Immunogen	Peptide sequence around site of 2-hydroxyisobutyryl-Lys (36) derived from Human Histone H2A type 1
Target	HIST1H2AG
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P0C0S8
Dilution Range	ICC:1:200-1:500
Expiration Date	12 months from date of receipt.



Immunocytochemistry analysis of Hela cel...