



Product Datasheet

Acetyl-HIST1H2AG (K15) antibody (orb517085)

Description	Acetyl-HIST1H2AG (K15) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF
Immunogen	Peptide sequence around site of Acetyl-Lys (15) derived from Human Histone H2A type 1
Target	HIST1H2AG
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P0C0S8
Dilution Range	ICC:1:20-1:200, IF:1:50-1:200
Expiration Date	12 months from date of receipt.



Immunocytochemistry analysis of HeLa cell...



Immunofluorescence staining of HeLa cell...