



Product Datasheet

HIST1H2AG (Ab-118) antibody (orb517097)


Description	HIST1H2AG (Ab-118) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, IP, WB
Immunogen	Peptide sequence around site of Lys (118) derived from Human Histone H2A type 1
Target	HIST1H2AG
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P0C0S8
Dilution Range	WB:1:500-1:5000, IHC-P:1:10-1:100, IF:1:1-1:10, IP:1:200-1:2000
Expiration Date	12 months from date of receipt.



Immunohistochemistry of paraffin-embedde...



Western blot aanalysis of NIH/3T3 whole ...



Immunoprecipitation analysis of NIH/3T3 ...