



Product Datasheet

MEK1/MEK2 (Phospho-Ser217/Ser221) Antibody (orb5171)

Description	MEK1/MEK2 (Phospho-Ser217/Ser221) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q02750
Expiration Date	12 months from date of receipt.



Western blot analysis
of UV
MEK1/MEK2(Ph...



Immunohistochemical
analysis of formalin...