



Product Datasheet

Acetyl-HIST1H4A (K31) antibody (orb517111)


Description	Acetyl-HIST1H4A (K31) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ChIP, ELISA, ICC, IF, WB
Immunogen	Peptide sequence around site of Acetyl-Lys (31) derived from Human Histone H4
Target	HIST1H4A
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P62805
Dilution Range	WB:1:500-1:5000, ICC:1:20-1:200, IF:1:10-1:100
Expiration Date	12 months from date of receipt.



Chromatin
Immunoprecipitation
Hela (4*10...



Western blot
analysis of Hela
whole cel...



Immunofluorescence staining of Hela cell...