



Product Datasheet

HIST1H1C (Ab-45) antibody (orb517131)


Description	HIST1H1C (Ab-45) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ChIP, ELISA, IHC, WB
Immunogen	Peptide sequence around site of Lys (45) derived from Human Histone H1.2
Target	HIST1H1C
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P16403
Dilution Range	WB:1:500-1:5000, IHC-P:1:10-1:100
Expiration Date	12 months from date of receipt.



Chromatin Immunoprecipitation
Hela (4*10...



Western blot aanalysis
of Hela whole cel...



Immunohistochemistry
of paraffin-embedde...