


Product Datasheet

Histone H3 (Ab-42) antibody (orb517144)

Description	Histone H3 (Ab-42) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Peptide sequence around site of Arg (42) derived from Human Histone H3.1
Target	HIST1H3A
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P68431
Dilution Range	IHC-P:1:10-1:100
Expiration Date	12 months from date of receipt.



Immunohistochemistry of paraffin-embedde...