



## Product Datasheet

### Histone H3 (Ab-128) antibody (orb517147)


|                            |                                                                           |
|----------------------------|---------------------------------------------------------------------------|
| <b>Description</b>         | Histone H3 (Ab-128) antibody                                              |
| <b>Species/Host</b>        | Rabbit                                                                    |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                         |
| <b>Conjugation</b>         | Unconjugated                                                              |
| <b>Tested Applications</b> | ELISA, IHC, IP, WB                                                        |
| <b>Immunogen</b>           | Peptide sequence around site of Arg (128) derived from Human Histone H3.1 |
| <b>Target</b>              | HIST1H3A                                                                  |
| <b>Preservatives</b>       | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol               |
| <b>Form/Appearance</b>     | Liquid                                                                    |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.             |
| <b>Note</b>                | For research use only                                                     |
| <b>Isotype</b>             | IgG                                                                       |
| <b>Clonality</b>           | Polyclonal                                                                |
| <b>Uniprot ID</b>          | <b>P68431</b>                                                             |
| <b>Dilution Range</b>      | WB:1:500-1:5000, IHC-P:1:10-1:100, IP:1:200-1:2000                        |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                           |



Immunohistochemistry of paraffin-embedde...



Western blot aanalysis of HeLa whole cel...



Immunoprecipitation analysis of HeLa who...