

Product Datasheet

DNAJC15 antibody (Biotin) (orb517344)

Description	DNAJC15 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human DnaJ homolog subfamily C member 15 protein (59-150AA)
Target	DNAJC15
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9Y5T4
Expiration Date	12 months from date of receipt.