

Product Datasheet

ATP2B4 antibody (Biotin) (orb517350)

Description	ATP2B4 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Plasma membrane calcium-transporting ATPase 4 protein (281-345AA)
Target	ATP2B4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P23634
Expiration Date	12 months from date of receipt.