

Product Datasheet

p63 antibody (FITC) (orb517358)

Description	p63 antibody (FITC)
Species/Host	Rabbit
Reactivity	Zebrafish
Conjugation	FITC
Immunogen	Recombinant Zebrafish Tumor protein p63 protein (41-153AA)
Target	tp63
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	E9QG65
Expiration Date	12 months from date of receipt.