

Product Datasheet

MGAT5 antibody (FITC) (orb517394)

Description	MGAT5 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Alpha-1,6-mannosylglycoprotein 6-beta-N-acetylglucosaminyltransferase A protein (107-217AA)
Target	MGAT5
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q09328
Expiration Date	12 months from date of receipt.