

Product Datasheet

PNCK antibody (FITC) (orb517448)

Description	PNCK antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Calcium/calmodulin-dependent protein kinase type 1B protein (1-115AA)
Target	PNCK
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q6P2M8
Expiration Date	12 months from date of receipt.