

Product Datasheet

CAND1 antibody (FITC) (orb517463)

Description	CAND1 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Cullin-associated NEDD8-dissociated protein 1 protein (1150-1230AA)
Target	CAND1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q86VP6
Expiration Date	12 months from date of receipt.