

## Product Datasheet

### MYO3B antibody (FITC) (orb517511)

|                        |                                                                  |
|------------------------|------------------------------------------------------------------|
| <b>Description</b>     | MYO3B antibody (FITC)                                            |
| <b>Species/Host</b>    | Rabbit                                                           |
| <b>Reactivity</b>      | Human                                                            |
| <b>Conjugation</b>     | FITC                                                             |
| <b>Immunogen</b>       | Recombinant Human Myosin-IIIb protein (1110-1276AA)              |
| <b>Target</b>          | MYO3B                                                            |
| <b>Preservatives</b>   | Preservative: 0.03% Proclin 300.<br>Constituents: 50% Glycerol   |
| <b>Form/Appearance</b> | Liquid                                                           |
| <b>Storage</b>         | Upon receipt, store at -20°C or -80°C.<br>Avoid repeated freeze. |
| <b>Note</b>            | For research use only                                            |
| <b>Isotype</b>         | IgG                                                              |
| <b>Clonality</b>       | Polyclonal                                                       |
| <b>Uniprot ID</b>      | <a href="#">Q8WXR4</a>                                           |
| <b>Expiration Date</b> | 12 months from date of receipt.                                  |