

Product Datasheet

WDR13 antibody (FITC) (orb517538)

Description	WDR13 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human WD repeat-containing protein 13 protein (14-156AA)
Target	WDR13
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9H1Z4
Expiration Date	12 months from date of receipt.