

Product Datasheet

CNOT10 antibody (FITC) (orb517595)

Description	CNOT10 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human CCR4-NOT transcription complex subunit 10 protein (595-744AA)
Target	CNOT10
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9H9A5
Expiration Date	12 months from date of receipt.