

Product Datasheet

2-hydroxyisobutyryl-HIST1H1C (K80) antibody (orb517631)

Description 2-hydroxyisobutyryl-HIST1H1C (K80) antibody

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications ChIP, ELISA

Immunogen Peptide sequence around site of 2-hydroxyisobutyryl-Lys (80) derived from Human Histone H1.2

Target HIST1H1C

Preservatives Preservative: 0.03% Proclin 300.
Constituents: 50% Glycerol

Form/Appearance Liquid

Storage Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.


Note For research use only

Isotype IgG

Clonality Polyclonal

Uniprot ID **P16403**

Expiration Date 12 months from date of receipt.



Chromatin
Immunoprecipitation
Hela (4*10...