



Product Datasheet

2-hydroxyisobutyryl-HIST1H1C (K158) antibody (orb517632)

Description	2-hydroxyisobutyryl-HIST1H1C (K158) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ChIP, ELISA, ICC, WB
Immunogen	Peptide sequence around site of 2-hydroxyisobutyryl-Lys (158) derived from Human Histone H1.2
Target	HIST1H1C
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P16403
Dilution Range	WB:1:500-1:5000, ICC:1:10-1:100
Expiration Date	12 months from date of receipt.



Chromatin
Immunoprecipitation
Hela (4*10...



Immunocytochemistry
analysis of Hela cel...