

Product Datasheet

ntpcr antibody (Biotin) (orb517679)

Description	ntpcr antibody (Biotin)
Species/Host	Rabbit
Reactivity	Zebrafish
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Zebrafish Nucleoside-triphosphatase, cancer-related protein (1-185AA)
Target	ntpcr
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	B8K0C4
Expiration Date	12 months from date of receipt.