

Product Datasheet

Ita antibody (HRP) (orb517681)

Description	Ita antibody (HRP)
Species/Host	Rabbit
Reactivity	Zebrafish
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Zebrafish Lymphotoxin alpha (TNF superfamily, member 1) protein (116-203AA)
Target	Ita
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q4W897
Expiration Date	12 months from date of receipt.