


Product Datasheet

rhoac antibody (orb517716)

Description	rhoac antibody
Species/Host	Rabbit
Reactivity	Zebrafish
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Zebrafish Rho-related GTP-binding protein RhoA-C protein (1-190AA)
Target	rhoac
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q7T2E8
Expiration Date	12 months from date of receipt.



Western blot
analysis of
Recombinant
pr...