


Product Datasheet

olig2 antibody (orb517720)

Description	olig2 antibody
Species/Host	Rabbit
Reactivity	Zebrafish
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Zebrafish Oligodendrocyte lineage transcription factor 2 protein (150-273AA)
Target	olig2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	A8WIK8
Expiration Date	12 months from date of receipt.



Western blot
analysis of
Recombinant
pr...