

Product Datasheet

CLEC14A antibody (HRP) (orb517773)

Description	CLEC14A antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human C-type lectin domain family 14 member A protein (278-349AA)
Target	CLEC14A
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q86T13
Expiration Date	12 months from date of receipt.