


Product Datasheet

TSPEAR antibody (orb517792)

Description	TSPEAR antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Thrombospondin-type laminin G domain and EAR repeat-containing protein (248-433AA)
Target	TSPEAR
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8WU66
Expiration Date	12 months from date of receipt.



Western blot analysis of Raji whole cell...