


Product Datasheet

GALNT9 antibody (orb517800)

Description	Rabbit polyclonal antibody to GALNT9
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Polypeptide N-acetylgalactosaminyltransferase 9 protein (184-338AA)
Target	GALNT9
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9HCQ5
Dilution Range	WB:1:500-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of SH-SY5Y whole ...