

Product Datasheet

ache antibody (FITC) (orb517806)

Description	ache antibody (FITC)
Species/Host	Rabbit
Reactivity	Zebrafish
Conjugation	FITC
Immunogen	Recombinant Zebrafish Acetylcholinesterase protein (37-155AA)
Target	ache
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9DDE3
Expiration Date	12 months from date of receipt.