


Product Datasheet

CCDC170 antibody (orb517840)

Description	CCDC170 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Coiled-coil domain-containing protein 170 protein (574-651AA)
Target	CCDC170
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8IYT3
Expiration Date	12 months from date of receipt.



Western blot
analysis of
HepG2 whole
ce...