

Product Datasheet

ARNTL antibody (Biotin) (orb517859)

Description	ARNTL antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Aryl hydrocarbon receptor nuclear translocator-like protein 1 protein (366-585AA)
Target	ARNTL
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O00327
Expiration Date	12 months from date of receipt.